



NEUROREADER®
BRINGS OBJECTIVITY INTO PRACTICE

Brainreader Mediation Server 3.3

Brainreader's mediation server installed on your site enables full integration with your facility's PACS. The mediation server 3.3 works with the newest Neuroreader® version: NR 2.7.



Features and Benefits

- ▶ Download The Neuroreader® Report as a PDF from the Mediation Server Web Interface.
- ▶ Multiple companies and accounts can use the same Mediation Server.
- ▶ Accessible link to send support information to Brainreader support.
- ▶ Password aging for better security.
- ▶ Ability to allow and deny filters (filters away incorrect DICOM files sent to the Mediation Server).
- ▶ New column displaying patient study ID for easier recognition of the MRI scan sent for analysis.
- ▶ Flexibility to send analysis reports as both encapsulated PDF and Secondary capture to PACS.
- ▶ NR 2.7 is customizable, which makes it possible to tailor the Neuroreader® Report to your exact needs:
 - ▶ The different sections can be freely moved to any desired location in the report.
 - ▶ Get more asymmetry indices (totalling 26)
 - ▶ Get 85 brain volumes (Left/Right of the Gray and White Matter volumes)
- ▶ Original and segmented 3D images in NIFTI and DICOM formats.
- ▶ 3D overlays in NIFTI and DICOM formats for:
 - ▶ Brain
 - ▶ Combined Classification
 - ▶ Filtered Combined Labels
 - ▶ Filtered Combined Labels in Color
 - ▶ Intracranial
- ▶ Exportable JSON values of the results for research purposes.

How it works

