



NEUROREADER®  
BRINGS OBJECTIVITY INTO PRACTICE

# Brainreader Mediation Server 3.3

Brainreader's mediation server installed on your site enables full integration with your facility's PACS. The mediation server 3.3 works with the newest Neuroreader® version: NR 2.7.




# Brainreader Mediation Server 3.3

## Features and Benefits

- ▶ Download The Neuroreader® Report as a PDF from the Mediation Server Web Interface.
- ▶ Multiple companies and accounts can use the same Mediation Server.
- ▶ Accessible link to send support information to Brainreader support.
- ▶ Password aging for better security.
- ▶ Ability to allow and deny filters (filters away incorrect DICOM files sent to the Mediation Server).
- ▶ New column displaying patient study ID for easier recognition of the MRI scan sent for analysis.
- ▶ Flexibility to send analysis reports as both encapsulated PDF and Secondary capture to PACS.
- ▶ NR 2.7 is customizable, which makes it possible to tailor the Neuroreader® Report to your exact needs:
  - ▶ The different sections can be freely moved to any desired location in the report.
  - ▶ Get more asymmetry indices (totalling 26)
  - ▶ Get 85 brain volumes (Left/Right of the Gray and White Matter volumes)
- ▶ Original and segmented 3D images in NIFTI and DICOM formats.
- ▶ 3D overlays in NifTI and DICOM formats for:
  - ▶ Brain
  - ▶ Combined Classification
  - ▶ Filtered Combined Labels
  - ▶ Filtered Combined Labels in Color
  - ▶ Intracranial
- ▶ Exportable JSON values of the results for research purposes.

## How it works



FR.29.8326.v1-US